

# Steps

Wednesday, 22 May 2024 4:37 PM

1. Wipe mastoids
2. Attach blue reference (left) and yellow ground (right), and attach to master SAHARA box
3. Put on cap
4. Measure front-back, center to cz
5. Measure left-right, center to cz
6. Look from front
7. Repeat measure to double check
8. Make sure all electrodes are standing up

=====

1. Power on and make sure dongle connected
2. Turn on Simulink model, two scopes, and check signals
  - a. Impedance check
  - b. Decrease bandpass to {0.5, 30} on one of the amps, check for ALPHA waves
  - c. Return bandpass to {0.5, 200}
3. Adjust electrodes